

Comparison of Various models of Hot Air Oven

Sr No.	PARTICULAR	HOT AIR OVEN		
	Model	HAO-225	HAO-95	HAO-27
1	Inner Chamber Dimensions (W X D X H)	24 " X 24 " X 24 " OR 61 X 61 X 61 cms.	18 " X 18 " X 18 " OR 46 X 46 X 46 cms.	12 " X 12 " X 12 " OR 30 X 30 X 30 cms.
2	Exterior Chamber Dimensions (A X B X C)	32 " X 32 " X 41" OR 81 X 81 X 104 cms.	26 " X 26 " X 35" OR 66 X 66 X 89 cms.	20 " X 20 " X 29" OR 50 X 50 X 73 cms.
3	Overall Dimensions (P X Q X R)			
4	Inner Chamber Capacity	225 Liters.	95 Liters.	27 Liters.
5	No Of Perforated Shelves	3 Nos (Adjustable at Different Height).		
6	Temperature Range	Ambient +10 °C to 200 °C.		
7	Temperature Resolution.	0.1 °C.		
8	Temperature Uniformity	± 6 % Of Set point Temperature.		
9	Input Supply Voltage	230 V,1 Phase,50 Hz.		
10	Air Circulation Fan	Provided with variable speed control.		